

Interventional
imaging

▶ SEAL ONE®

All-In-One Radial Compression Device

SEAL ONE



- Precise adjustment and control of compression level
- Targeting the radial artery
- Enhanced visibility of the puncture site
- Secured decompression

NEW

PEROUSE
MEDICAL

▶ SEAL ONE®



○ FEATURES

▶ Precise adjustment and control of compression level⁽¹⁾

Compression knob with scale

- Limited risk of radial artery occlusion
- Optimization of use time

▶ Targeting the radial artery

Specific design of the compression pad

- No compression of the ulnar artery
- Ensured venous return

▶ Enhanced visibility of the puncture site

Transparent device

- Better visual control

▶ Two-step decompression

Secured button

- Limited risk of unintended decompression

▶ Accessory free

Compression / Decompression knob

- Simplified handling

▶ Display of positioning time

Time indicator

- Easier monitoring of decompression protocol

| Reference | Description | Packaging |
|-----------|---------------------------|-------------|
| 0259NA | Radial compression device | 30 unit-box |

(1) Each half degree of graduation increases the compression by 0.35mm and reciprocally for the decompression.

