

Interventional
imaging

▶ MYSHELL LITE

The Y connector
designed for your hand



- Push-pull closing system
- Double silicone valve
- Inner diameter up to 9F



○ CONNECTION LINE AND STOPCOCK

- connection line
 - ▶ no displacement of the catheters
- stopcock
 - ▶ drug injection

○ BODY AND BUTTON

- push-push system
 - ▶ one hand used
- compact system
 - ▶ ergonomic handling

○ DOUBLE SILICONE VALVE

- wire and catheter control through closed valve
 - ▶ minimum handling
 - ▶ limited blood leaking and pressure drop
- soft material
 - ▶ keep the integrity of the coated catheters
- 7F / 9F option
 - ▶ suitable for kissing,...

| References Inner diameter 7F | References Inner diameter 9F | Description | Packaging |
|---------------------------------|---------------------------------|--|-------------|
| 0210PC | | myshell lite Y haemostatic valve 7F 20 cm PVC connection line - 3-way stopcock | 25-unit box |
| 0210TL | | PTCA set 3 items myshell lite 7F 20 cm PVC connection line - 3-way stopcock metallic guidewire introducer screwing torque device | 25-unit box |
| | 0210QA | myshell lite Y haemostatic valve 9F 20 cm PVC connection line - female Luer lock | 25-unit box |
| | 0210QC | myshell lite Y haemostatic valve 9F 20 cm PVC connection line - 3-way stopcock | 25-unit box |
| | 0210QG | myshell lite Y haemostatic valve 9F 50 cm PVC connection line - 3-way stopcock | 25-unit box |
| | 0210SE | PTCA set 3 items myshell lite 9F 20 cm PVC connection line - 3-way stopcock metallic guidewire introducer alligatork torque device | 25-unit box |

Please contact PEROUSE MEDICAL for non sterile and OEM products

