

Interventional
Imaging


▶ SCREWING VALVE

SCREWING

- Large screwing ring
- Opened axial angle
- Inner diameter up to 9F



▶ SCREWING VALVE



● SCREWING RING

- Large diameter
- ▶ easy screwing/unscrewing

● SILICONE VALVE

- Inner diameter 9F
- ▶ suitable for Rotablator®

● AXIAL LINE

- Opened angle
- ▶ ergonomic handling

Product code	Description	Packaging
0210RM	<ul style="list-style-type: none"> • Screwing valve 9F • Metallic guidewire introducer • Screwing torque device 	50-unit box
0218TR	<ul style="list-style-type: none"> • Flamingo 30 atm inflation device • PTCA set 3 items • Screwing valve • Metallic guidewire introducer • Screwing torque device 	12-unit box
0185TS	<ul style="list-style-type: none"> • Dolphin 30 atm inflation device • PTCA set 3 items • Screwing valve • Metallic guidewire introducer • Screwing torque device 	6-unit box

